

KANGAROO™ and ARGYLE™ neonatal and pediatric feeding tube

Features and Benefits

- Purple ENFit™* feeding port indicates nutritional port for enhanced safety and reduced risk of wrong route errors.
- DEHP-free PVC or Thermosensitive polyurethane material options available
- Depth markings at 1cm intervals facilitates tracking of tube migration to improve patient safety
- Radiopaque, purple line
- Sterile for optimal patient protection
- Size available 3.5FR, 5FR, 6.5FR, 8FR, 10FR
- Length available 31cm, 51cm, 91cm, 107cm, 41 cm (for Argyle™ PVC only)

The new standard in enteral feeding connections

Intravenous and enteral systems previously had very similar connectors. A misconnection could result in enteral formula being delivered into a vein or intravenous medication being delivered into the stomach. These misconnection occurrences, can cause serious health risks.

In order to avoid misconnections, an international group of clinicians, manufacturers and regulators have developed a new standard ISO 80369-3. The objective is to improve patient safety by reducing the incidence of medical device tubing misconnections. As a result all new enteral medical devices will carry the new ENFit™* connection design.



Connecting for greater patient safety



GEDSA

Cardinal Health is a chartered member of GEDSA. The Global Enteral Device Supplier Association (GEDSA) was formed to help introduce international standards in medical device tubing connectors, intended to enhance patient safety. www.gedsa.org

Stay connected

Stay Connected is a communications program created by GEDSA to facilitate a successful transition to new, ENFit™* connectors intended to promote compatibility and consistency while reducing the likelihood of tubing misconnections. The Stay Connected initiative will help gradually introduce new standard connectors for specific delivery systems, including neuraxial, limb cuff inflation, and respiratory applications.

Sign up to stay connected

To sign up for email updates with the latest information and tools to help you with this transition, visit: www.stayconnected.org †

Code	Description	Case QTY
------	-------------	----------

Kangaroo™ Polyurethane Feeding Tubes

461355E	3.5 Fr – length 31 cm	10
461412E	5 Fr – length 51 cm	10
461420E	6.5 Fr – length 51 cm	10
461438E	6.5 Fr – length 91 cm	10
461701E	8 Fr – length 51 cm	10
461503E	5 Fr – length 91 cm	10
461800E	8 Fr – length 107 cm	10
461042E	10 Fr – length 107 cm	10

Argyle™ PVC Feeding Tubes

461206E	3.5 Fr – length 31 cm	10
460802E	5 Fr – length 41 cm	50
461008E	5 Fr – length 91 cm	50
460711E	6.5 Fr – length 41 cm	50
460604E	8 Fr – length 41 cm	50
460406E	8 Fr – length 107 cm	50
460208E	10 Fr – length 107 cm	50

† This document may provide links or references to third party websites. Cardinal Health has not reviewed any third party website and is not responsible for the content or the privacy policies of any third party website and shall not be liable for any damages or injury arising from that content. Any links to other sites are provided as merely a convenience.